



The Comfort



The Handling of heavy obese patients is a serious health and safety matter. Adopting the best equipment and practices will ensure comfort and security. The Bariatric Conveyance Chair has been designed to meet the expressed needs of clinical providers. The Chair offers full support to all upper body soft tissue as well as being wide enough to give full base comfort. All the surfaces are antimicrobial with 50mm foam cushions and over the arm rests. The seat height is 500mm to allow easy access for short as well as tall patients. The chair can be reclined to 45deg to give the patient a sense of added security during transportation and forward 10deg. to assist in seating and standing.

Features

- Multi-functional treatment and conveyance chair for leg and head procedures
- Reclines to 45 degrees and tips forward for easy patient seating and standing
- Deep cushioned seat and back with anti-microbial upholstery
- Motor driven leg rest with treatment platforms
- Colour range of 54 options
- Duplex castors for easy movement
- Optional 5th wheel electric drive

Specifications

- Easy to move transport chair
- Large adjustable footrest with separate treatment surfaces
- Reclines to 45deg. and tips forward 15 deg.
- Rigid sturdy and cushioned armrests for support
- 50mm deep upholstery with profiled surface
- Seat widths of 940mm and 700mm
- Tip forward of 10 deg. to help seating and standing
- Low seat height to suit all patients
- One lockable wheel for easy steering
- Centre mounted drive wheel to give easy operation as an option
- Centre drive wheel can be lowered and raised
- Integral battery charger
- Fingertip control when the drive wheel is operated
- SWL 1000lbs (450kg)
- Battery operated with integral charger
- Sturdy welded steel frame
- Complies with the Medical Directive 93/94/EEC



Installation

- No installation or building works required

sales@designspecific.co.uk

Design Specific Ltd.
Unit 10
Caburn Enterprise Park
The Broyle
Ringmer
BN8 5NP
United Kingdom

T: +44(0)1273 813904
M: +44(0)7501 464703